

FOX

Y-CONNECTOR



Features

- click closing system
- three-stage valve creates maneuverability without any bleeding
- auto-closing mechanism ensures the device to move from semi-closed position to closed position automatically when contrast is injected from the side-on lumen
- ergonomic design allows for one-hand operation

Product details

- inner diameter up to 8F
- three-stage valve
- pressure line (20 cm) and stopcock
- ergonomic handling

REF. NUMBER	DESCRIPTION
DRG-HV-8CP-TB	Y-connector click type with pressure line and stopcock