

FALCON

INFLATION DEVICE



Features

- ergonomic design facilitates doctor's operation
- trigger release mechanism provides easy and quick deflation
- Inclination of manometer provides easy and precise view of pressure points

Product details

- barrel volume: 25ml
- pressure range: 30atm/bar
- 3-way stopcock, sterile, tyvek-plastic blister
- spiral design ensures precise pressure
- overpressure protective device provides safer pressure control
- release handle makes instant and easy inflation and deflation

| REF. NUMBER | DESCRIPTION |
|------------------|---|
| DRG-ID-A2530H-TB | 25ml, 30atm/bar, Stopcock, Sterile, Tyvek-plastic Blister |